# Aquasil<sup>®</sup> LV

## Smart Wetting<sup>®</sup> Impression Material

### DIRECTIONS FOR USE - ENGLISH

**Caution: For dental use only.**  
U.S.A.: Rx only.

#### 1. PRODUCT DESCRIPTION

Aquasil LV Smart Wetting<sup>®</sup> Impression Material is a quadrifunctional hydrophilic addition reaction silicone, light body, elastomeric impression material with excellent hydrophilic properties, dimensional accuracy, high tear strength, and resistance to permanent deformation. Aquasil LV Impression Material with Smart Wetting<sup>®</sup> improves both wetting on the tooth surface and model detail reproduction.

#### 1.1 Delivery forms\*

Aquasil LV Impression Material is available in:

- 50mL cartridges for dispensing gun
- 90mL Tubes for manual mixing

\*Some delivery forms may not be available in all countries.

#### 1.2 Composition

Polydimethylsiloxane polymer; Polymethylhydrogen Siloxane; Calcium Sulfate; Silicon Dioxide; Sodium Aluminosilicate; Organic Platinum Complex; Surfactant; Titanium Dioxide; Metallic Oxide Pigments

#### 1.3 Indications for use

Aquasil LV Impression Material is suitable for:

- All impression techniques where a light body (wash) material would be desired by the operator.
- Dual phase, one-step ("double-mix") impression technique with Aquasil Tray Impression Materials.
- Dual phase, two-step ("putty-wash") impression technique with Aquasil Putty (all versions) or Reprosil<sup>®</sup> Putty Impression Materials.
- Aquasil LV Fast Set Impression Material is designed to be used with Aquasil Rigid Fast Set Impression Material in a dual phase Impression Technique, capturing one preparation only (single unit crown).
- Dual phase, for the duplication of models.

#### 1.4 Contraindications

None known.

#### 1.5 Technical specifications

Classification: ISO 4823, Type 3, Light-Bodyed Consistency

Working Time (22°C/72°F)..... **Fast Set:** 1'15" - 1'45"  
**Regular Set:** 2'15" - 2'45"  
(room temperature mixing and preparation, maximum insertion time)

Working Time (37°C/98°F)..... **Fast Set:** 35" Maximum  
**Regular Set:** 1'10" Maximum  
(intraoral syringing prior to tray insertion)

Setting Time (Mouth Removal Time)..... **Fast Set:** 3'00" from start of mix  
**Regular Set:** 6'00" from start of mix

% Linear Dimensional Change..... < 0.50  
% Recovery from deformation..... ≥ 98  
Detail Reproduction..... .20 microns

#### 1.6 Compatible tray adhesive

Aquasil Impression Material is compatible with conventional tray adhesives such as Caulk<sup>®</sup> Tray Adhesive or Silfix<sup>®</sup> Tray Adhesive designed for use with vinyl polysiloxane materials (See complete Directions for Use).

## 2. GENERAL SAFETY NOTES

Be aware of the following general safety notes and the special safety notes in other chapters of these directions.

#### Safety alert symbol

This is the safety alert symbol. It is used to alert you to potential personal injury hazards. Obey all safety messages that follow this symbol to avoid possible injury.

#### 2.1 Warnings

Aquasil Impression Material contains polymerizable monomers and surfactants which may be irritating to skin, eyes and oral mucosa and may cause allergic contact dermatitis in susceptible persons.

- **Avoid Eye contact** to prevent irritation and possible corneal damage. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention.
- **Avoid Skin contact** to prevent irritation and possible allergic response. In case of contact, reddish rashes may be seen on the skin. If contact with skin occurs immediately remove material with cotton and wash thoroughly with soap and water. In case of skin sensitization or rash, discontinue use and seek medical attention.
- **Avoid ingestion/swallowing** of material to prevent irritation or obstruction. If accidental swallowing occurs, drink lots of water. This material is not hazardous when small quantities are ingested. Larger quantities may cause bowel obstruction. Seek medical attention in the event of digestive irregularities.
- **Do not use Aquasil Impression Material as a temporary reliner.** Use of the material in this manner may cause irritation to the oral mucosa. If symptoms occur discontinue use immediately. If irritation does not subside consult a physician.
- **Avoid multiple impressions at one sitting** (more than three of the same arch/area) to prevent irritation to the oral mucosa. If symptoms occur discontinue use immediately. If irritation does not subside consult a physician.

#### 2.2 Precautions

- This product is intended to be used only as specifically outlined in the Directions For Use. Any use of this product inconsistent with the Directions For Use is at the discretion and sole responsibility of the practitioner.
- Insufficient data exist to support the use of Aquasil Impression Material for impressions on exposed bone. Aquasil Impression Material should only be used on or around intact mucosa.

- Devices marked "single use" on the labeling are intended for single use only. Discard after use. Do not reuse in other patients in order to prevent cross-contamination.
- All products should be used at room temperature. Higher temperatures reduce work times and laboratory bench set times (faster), lower temperatures increase them (slower).
- To obtain optimum physical properties, the mix should be completely homogeneous (streak free). Variations in the relative amounts of base and catalyst will neither significantly increase nor decrease work or set times. Optimal performance requires a ratio of 1:1.
- Aquasil Impression Material should be used with properly fitted, rigid trays for most applications. Trays with retentive features such as Caulk Rim-Lock<sup>®</sup> Trays or perforations, pretreated with tray adhesive, are recommended. To reduce gagging, choking or swallowing, do not overfill tray.
- Material should extrude easily: DO NOT USE EXCESSIVE FORCE. Excessive pressure may result in unanticipated extrusion of the material or cause cartridge rupture.
- Aquasil LV Fast Set is not recommended for full arch impression or multiple units. Syringing more than two preparations may cause operator to exceed the allowable work time, causing incomplete and/or inaccurate impression.

#### Interactions:

- Do not use Hydrogen peroxide-containing dental materials in conjunction with this product since they may interfere with the hardening of the polymeric components.
- Use of latex or "sulfur-containing" polymeric gloves, some hand soaps and lotions can lead to interference with setting reaction. Contamination may accumulate on the surface of the skin and may not easily wash off. Non-latex gloves worn during use will avoid this possibility.
- Handle cords with college pliers. Do not touch prepared tooth with gloves. Rinse preparation with copious amounts of water before making the impression.
- Aquasil Impression Material should not be intermixed with, or used in conjunction with, other manufacturer's vinyl polysiloxanes or any polyether, polysulfide or conventional (condensation cured) silicones.

#### 2.3 Adverse reactions

- **Eye contact:** irritation and possible corneal damage.
- **Skin contact:** irritation or possible allergic response. Reddish rashes may be seen on the skin.
- **Bowel obstruction** or other digestive distress may result from ingestion of mixed impression material. (See Warnings)

#### 2.4 Storage conditions

Inadequate storage conditions may shorten the shelf life and may lead to malfunction of the product.

- Store at temperatures between 10°C/50°F - 24°C/75°F.
- Use the product at room temperature.
- Do not store used cartridges with new (unused) mixing tips attached.
- Close tubes immediately after use. Do not interchange base and catalyst screw caps.
- Do not freeze.
- Protect from moisture.
- Do not use after expiration date.

## 3. STEP BY STEP INSTRUCTIONS

### 3.1 Cartridge Dispensing Gun

1. Insert cartridge into dispensing gun.
2. Bleed cartridge before installing new mixing tip.
3. Preparation of syringe/wash material
  - Dispense directly to the mouth with the intraoral tip; **OR**
  - An intraoral syringe may be filled with mixed material from the cartridge.
4. Disinfect cartridge dispensing gun, cartridge and attached used mixing tip for storage.
  - Immediately prior to subsequent usage, remove the used mixing tip, bleed cartridge, install a new mixing tip.
  - Dispose of cartridge when empty. Do not attempt to refill or consolidate used cartridges.

### 3.2 Tube Dispensed

1. Dispense equal lengths of Aquasil LV Impression Material catalyst and base onto mixing pad.
2. Using a stiff spatula and a stropping action, mix impression material. To obtain optimum physical properties, the mix should be completely homogeneous (streak-free). Spatulate for approximately 30-45 seconds mix time.
3. Backfill the impression syringe and/or load impression tray.

### 3.3 Impression technique: Dual Phase one-step ("double-mix") Technique with a Rigid Tray or with Double Arch Technique

**NOTE: Tray and wash material mix timing:** It is highly recommended to utilize an assistant to assure simultaneous syringing and tray loading.

1. Assure hemostasis and clean the field with air water spray.
2. Syringe material: Backfill the impression syringe or attach an intraoral tip to the end of the mixing tip.  
Tray material: Dispense material (assistant) into impression tray. Do not overfill tray.
3. Simultaneous with tray loading, syringe (dentist) completely around the preparation/abutment.  
**NOTE:** Tray and wash material insertion timing
  - **Fast Set** wash and tray material, loaded tray must be inserted within 35 seconds from time of first syringing of wash material intraorally.
  - **Regular Set** wash and tray material, loaded tray must be inserted within 1 minute 10 seconds from time of first syringing of wash material intraorally.
4. Insert filled tray into the mouth.
5. The Minimum Removal Time is 3 minutes (Fast Set), 6 minutes (Regular Set) from the start of mix.
6. Follow disinfection procedures before shipping and/or pouring.

#### Timing Guide

	<b>Fast Set:</b> Up to 0'35" (35 sec) <b>Regular Set:</b> Up to 1'10" (1 min 10 sec)	<b>Fast Set:</b> 3'00" (3 min) <b>Regular Set:</b> 6'00" (6 min)
0:00 (0 sec)		
Begin simultaneous syringing and tray loading	Complete syringing and seat loaded tray	Remove set impression

**3.4 Impression technique:** Dual Phase two-step (“putty-wash”) Technique: With a primary tray impression.

1. Complete primary impression with Aquasil Putty or Reprosil® Putty material.
2. It is recommended to utilize an assistant to assure simultaneous syringing and tray loading.
3. Follow impression steps above, using wash material for both intraoral syringing and loading into primary tray impression.

**3.5. Impression Technique:** Duplication of Models

1. Load the impression tray with Aquasil Monophase or Aquasil Rigid tray material.
2. It is suggested to simultaneously syringe, using Aquasil LV Impression Material, completely around the model teeth.
3. Seat model into the loaded tray. Retain impression in position until firmly set.
4. Remove model from impression by pulling slowly to break seal.

#### 4. HYGIENE

##### Cross-contamination

- Do not reuse single use products. Dispose in accordance with local regulations.
- Reprocess reusable products as described below.

##### 4.1 Cleaning the dispenser

The cartridge dispensing guns may be cleaned by scrubbing with hot water and soap or detergent. Excess material may be cleaned with alcohol-moistened gauze.

##### 4.2 Disinfection of the dispenser gun and cartridge and tubes

Disinfect dispenser guns, cartridges and used tips left in place on cartridges for storage and Tubes with a hospital-level tuberculocidal disinfectant according to national and local regulations. Use impregnated wipes not sprays. **NOTE:** Destruction of the label by vigorous wiping. Wipe cartridge gently.

Some phenolic based agents and iodophor based products may cause surface staining.

Agents containing organic solvents, such as alcohol may tend to dissolve the plastic dispenser and plunger. The disinfectant manufacturer’s directions should be followed properly for optimum results.

**Technique Tip:** Unset Aquasil Impression Material can be removed from clothing with dry cleaning solvent.

##### 4.3 Disinfection of the Impression Material

The impression should be disinfected with a hospital-level tuberculocidal disinfectant according to national and local regulations.

- Iodophors, sodium hypochlorite (0.525%), chlorine dioxide, and dual or synergized quaternary ammoniums are approved disinfectants.
- Thoroughly soak by spraying or immerse the impression in any recommended hospital level disinfectant for the contact time recommended by the disinfectant manufacturer for optimum results.
- Water-based disinfectant solutions are preferred.

##### 4.4 Casting

- The impression should be removed from the disinfectant, rinsed with water and lightly air dried.
- The impression may be poured in plaster or stone 15 minutes after disinfection or delayed up to fourteen days. For epoxy dies, delay pour 60 minutes.
- Do not store impression in direct sunlight.
- If the impression is to be shipped, use suitable packaging to preclude distortion.
- The material is compatible with a range of dental plasters and die stones. Impressions may be silver or copper-plated.

#### 5. LOT NUMBER AND EXPIRATION DATE

1. Do not use after expiration date. ISO standard uses: “YYYY/MM.”
2. The following numbers should be quoted in all correspondences:
  - Reorder number
  - Lot number
  - Expiration date

Manufactured by  
Dentsply Caulk  
38 West Clarke Avenue  
Milford, DE 19963 USA  
Made in USA  
1-302-422-4511  
www.dentsplysirona.com

Dentsply DeTrey GmbH  
De-Trey-Strasse 1  
78467 Konstanz  
Germany  
49-7531-583-0  
www.dentsplysirona.com

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