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# Caulk® 34% Tooth Conditioner Gel Phosphoric Acid Gel

### DIRECTIONS FOR USE – ENGLISH

**CAUTION:** U.S. FEDERAL LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A DENTIST.

**Caulk® 34% Tooth Conditioner Gel** is a surface-conditioning agent for treatment of enamel and/or dentin prior to adhesive restoration placement. The product is supplied in syringes and bottles.

**COMPOSITION**  
Water, 34% Phosphoric acid, Silicon Dioxide, Surfactants, Blue Colorant

**INDICATION**  
Caulk® 34% Tooth Conditioning Gel (syringe or bottle dispensed) is used to surface condition enamel (and dentin where the total etch technique is desired) prior to restoration with adhesive restorative materials.

**CONTRAINDICATION**  
Caulk® 34% Tooth Conditioning Gel is contraindicated for use on dental pulp tissue.

	<b>Safety alert symbol</b>
<b>CAUTION</b>	This is the safety alert symbol. It is used to alert you to potential personal injury hazards. Obey all safety messages that follow this symbol to avoid possible injury.

### WARNINGS

- Caulk® 34% Tooth Conditioning Gel contains phosphoric acid, which may cause burns. Avoid contact with oral tissues, eyes and skin. Do not take internally.
- Eye Contact:** Caulk® 34% Tooth Conditioning Gel contains phosphoric acid, which may cause burn or may be irritating to eyes. Before using this product wear protective glasses as well as covering the patient's eyes to protect from squirting the material inadvertently. In case of accidental contact with eyes, rinse eyes immediately with plenty of water and seek medical attention.
- Skin Contact:** Caulk® 34% Tooth Conditioning Gel contains phosphoric acid, which can cause burns or skin sensitization in susceptible individuals. If contact with skin occurs immediately wipe off and flush with generous amounts of water, then wash well with soap and water after contact. If skin rash, irritation, sensitization or other allergic reaction occurs, seek medical attention immediately.
- Oral Mucosa Contact:** Avoid contact with oral soft tissues. If accidental contact occurs, flush mucosa with plenty of water and exporatorate water. If sensitization of mucosa persists, seek medical attention immediately.
- Ingestion:** Do not swallow or take internally. If accidental swallowing occurs, drink lots of water. If nausea or illness develop, seek medical attention immediately. Contact regional Poison Control Center immediately.
- Caulk® 34% Tooth Conditioning Gel in syringes should extrude easily: DO NOT USE EXCESSIVE FORCE. Excessive pressure may result in unanticipated extrusion of the material.
- For full coverage vital crown preparations, condition remaining enamel only. Etching of full coverage dentin surfaces is not recommended to minimize the possibility of post-operative sensitivity.
- Caulk 34% Tooth Conditioning Gel should not be used with patients who have a history of severe allergic reaction to any of the components.

### PRECAUTIONS

- This product is intended to be used only as specifically outlined in the Directions For Use. Any use of this product inconsistent with the Directions For Use is at the discretion and sole responsibility of the practitioner.
- Wear suitable protective eyewear, clothing and gloves. Protective eyewear is recommended for patients.
- Devices marked "single use" on the labeling are intended for single use only. Discard after use. Do not reuse in other patients in order to prevent cross-contamination.
- Replace original cap of Caulk® 34% Tooth Conditioning Gel tightly after each use to avoid evaporation. Discard dispensing tip after use, as dispensing tips may clog if gel is allowed to dry inside.
- DENTSPLY Caulk supplies the appropriate dispensing tip for the Caulk® 34% Tooth Conditioning Gel. This tip is the only tip that should be utilized for placement of the gel.
- For use on full coverage vital crown preparations see Warnings and Step-by-Step Instructions.
- The efficacy of Caulk® 34% Tooth Conditioning Gel as an etchant/cleanser of indirect restorations has not been demonstrated.
- Storage and Shelf Life:** Store syringe with original cap only. Do not store syringe with dispensing tip in place, as this may harden the material making extrusion difficult (see Warnings). Not to be stored at temperature exceeding 25°C. Do not freeze. Do not use after expiration date.
- Use only in well ventilated areas.

### ADVERSE REACTIONS

- Product may irritate the eyes and skin. **Eye contact:** irritation and possible corneal damage. **Skin contact:** irritation or possible allergic response. Reddish rashes may be seen on the skin. **Mucous Membranes:** inflammation, edema, sloughing. See Warnings.
- Product may cause serious health effects if ingested. See Warnings.
- Product may cause pulpal effects. See Contraindication and Warnings.
- Inhalation of vapors may cause varying degrees of damage to the affected tissue and also increased susceptibility to respiratory illness. See Precautions.
- The following medical conditions are generally aggravated by exposure: Individuals with pre-existing skin disorders, eye problems or impaired liver, kidney and respiratory functions.

### STEP-BY-STEP INSTRUCTIONS FOR USE I. CLEANING

Clean uninstrumented enamel with a rubber cup and water/flour of pumice or a cleaning paste such as Zircate®. Wash thoroughly with clean, oil-free water spray.

### II. PULP PROTECTION

In deep cavities, cover the dentin close to the pulp (less than 1 mm remaining) with a hard setting calcium hydroxide liner (Dycal® Radiopaque Calcium Hydroxide Composition) leaving the rest of the cavity surface free for bonding.

### III. APPLICATION

**Dispensing From Syringe:**

- Remove cap from the end of the syringe. To assure free flow of material from syringe, express a small amount onto a mixing pad.
- Attach disposable dispensing tip to the end of the syringe. Turn tip clockwise until fully seated. Tug on the tip to assure locking into syringe collar.
- Dispensing tip may be bent for easy access. Gentle curving of the tip is preferable to an acute bend. Bending may occlude the tip. Be sure to check that this has not occurred by extruding a small amount through the tip.

- Gel should flow freely with GENTLE PRESSURE. DO NOT USE EXCESSIVE FORCE. If not, remove syringe from patient field and check for obstruction.
- Apply gel through syringe dispensing tip to surfaces to be conditioned.

#### Dispensing From Bottle:

- Dispense desired amount onto a mixing pad or dappen dish. To prevent the gel from drying out, dispense gel just prior to placement procedures. Stirring will reduce viscosity.
- Using a clean, soft brush, apply gel to surfaces to be conditioned.

#### IV. CONDITIONING TOOTH SURFACE

- Apply Caulk' 34% Tooth Conditioning Gel to clean enamel. Allow conditioner to remain in place for 15 seconds. Where total etch technique is desired, begin with application of conditioner to enamel first, then flowing onto dentin. Allow conditioner to remain in place for a total of 15 seconds from initial enamel contact.
- Rinse conditioned areas thoroughly with water (approx. 10-15 seconds).
- If only enamel has been conditioned, air dry with clean, oil-water-free compressed air for 15 seconds. Properly conditioned enamel will appear dull, frosty white. **NOTE:** Protect exposed dentin from air to prevent desiccation.
- If total etch technique has been employed, blot dry conditioned areas with a moist cotton pellet. Do not rub tooth surface while drying. Properly conditioned dentin should appear glistening, without pooling of moisture.
- Proceed immediately to bonding technique desired.

Properly conditioned surfaces must be kept free of contamination by oral fluids. If salivary contamination occurs, thoroughly clean with forceful water spray (free of oil contamination), blot dry and repeat conditioning procedure on enamel only for 10 seconds, rinse as in step 2 above and dry as in step 3 above.

#### HYGIENE

#### Cross-contamination

- Do not reuse single use products. Dispose in accordance with local regulations.
- Reprocess reusable products as described below.

#### DISINFECTION OF SYRINGE AND/OR BOTTLE

The syringe and/or bottle if exposed spatter or spray of body fluids or that may have been touched by contaminated hands, or oral tissues, should be disinfected with a hospital-level disinfectant. Acceptable disinfectants are EPA-registered as tuberculocidal. Iodophors, sodium hypochlorite (5.25%), chlorine dioxide and dual or synergized quaternaries are approved disinfectants. Disinfect the syringe and/or bottle by spraying the syringe and/or bottle with any recommended hospital-level disinfectant except neutral glutaraldehyde for the contact time recommended by the disinfectant manufacturer for optimum results. Some phenolic-based agents and iodophor-based products may cause surface staining. Agents containing organic solvents, such as alcohol, may tend to dissolve the plastic. The disinfectant manufacturer's directions should be followed properly for optimum results. Water-based disinfectant solutions are preferred.

#### LOT NUMBER AND EXPIRATION DATE

- Do not use after expiration date. ISO standard uses: "YYYY/MM."
- The following numbers should be quoted in all correspondences.
  - Reorder number
  - Lot number on the syringe or bottle
  - Expiration date

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All rights reserved.
Dentsply DeTrey GmbH