SmartCem®2 Selbstadhäsiver Zement

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3.3 Cimentation technique			
	Self-Cure	Dual Cure – Non-Light Transmissible Restorations	Light Cure – Light Transmissible Restorations
TYPE OF RESTORATION	all	PFM, Zirconia, Alumina, opaque ceramics & composites	Translucent ceramics & composites
1. FILL AND SEAT CROWN	up to 2 min	up to 2 min	up to 2 min
<i>Work time extraoral from start of mix:</i>			
<i>Protect restoration from movement during gel phase cleanup through final set</i>			
2. REMOVE EXCESS			
Gel stage from seating:	1-2 min	up to 10 sec light cure	up to 10 sec light cure
Gel Duration:	45 sec		45 sec
3. STABILIZE			
Constant occlusal pressure	Constant	Light cure margins 20-40 sec	Light cure margins 20-40 sec.
4. FINAL SET			
6 min from start of mix	6 min from start of mix		Upon completion of light curing each surface" 10 sec <i>"Buccal, Lingual and Occlusal</i>

Danger of Injury due to excessive force

- Apply slow and steady pressure on the syringe
- Do not use excessive force – dual barrel syringe rupture may result

3.3.1 Dual Barreled syringe dispensed

1. Remove syringe cap. Dispense and discard a small amount of material from the dual-barreled syringe. Be sure material is flowing freely from both ports. Holding syringe vertically, carefully wipe away excess from base and catalyst do not cross contaminate and cause obstruction of the ports. Save syringe cap for replacement following use.

2. Install a mixing tip on the cartridge by lining up the v-shaped notch on the outside of the mixing tip with the v-shape notch on the syringe flange. Turn colored mixing tip cap 90 degrees in a clockwise direction to lock in place on syringe.

3. Gently depress syringe plungers to begin the flow of material. DO NOT use EXCESSIVE FORCE. If force is encountered, remove syringe from operating field, remove and discard mixing tip. Check for obstruction and confirm material is flowing from both syringe barrels. Wipe barrels and install new mixing tip as outlined above. Dispense a small amount through the mixing tip onto a mixing pad and discard.

4. Without delay, using gentle pressure, apply a thin, uniform layer of cement to the entire internal surface of the restoration directly from the mixing tip. At room temperature, SmartCem™ 2 Cement offers a minimum work time of 2 minutes. **Technique Tip:** The mixing tip may be bent slightly to allow direct intraoral access for placement of cement into preparations with internal anatomy. For endodontic post spaces, use of a Lentulo Spiral or metal file to aid placement in the post space is recommended.

5. Immediately seat the restoration in the mouth. Verify compete seating. A gentle rocking or vibratory motion may be helpful to insure optimal seating.

6. Protect restoration from contamination and movement until the final set of the cement (6 minutes from start of mix or in the case of light-transmissible restorations, completion of light curing).

3.4. Cleaning marginal excess SmartCem™ 2 Cement

3.4.1 Self-cure cleanup

The excess cement will reach the “gelled” state after approximately 1-2 minutes in the mouth, allowing easy removal. Excess cement will remain in the “gelled” state for approximately 1 minute. If exposed to directed operating light, “gel” state may be reached sooner and remain “gelled” for a shorter period. Immediately after reaching the “gelled” state, floss interproximally, only in the direction of the tooth seating, to remove excess cement. Complete excess cement removal using an instrument such as a rubber tip, a scaler or an explorer. **NOTE:** Cement within the crown has not yet set. Do not move, torque, or disturb the crown during cleanup.

Technique Tip: Following all excess removal, exposed margins may be light cured 20-40 seconds to assist restoration stabilization.

3.4.2 Optional dual-cure cleanup

Due to the dual-cure property of SmartCem™ 2 Cement, the operator has the option of utilizing a curing light to facilitate cleanup. Light curing to facilitate cleanup must be accomplished within the first minute following intraoral insertion. Excess cement cleanup may begin immediately following a brief exposure with the curing light. Conventional quartz tungsten halogen or LED lights producing light of wavelength 470nm are recommended. Use of narrow spectrum output lights or high power output may produce unexpected results. *Light curing mixed cement after one minute, or continuously for more than 10 seconds, at any time, will cause cement to adhesively set, making cleanup difficult.*

Immediately after seating is verified, briefly light-cure excess cement at the margins by constantly moving the curing light up around the margins for no more than 10 seconds. Excess cement will reach a “gelled” state after this 10 seconds cure.

The excess cement will remain in the “gelled” state for approximately 45 seconds following light exposure. All excess cement must be removed before final self-cure set is achieved, as outlined above.

NOTE: Cement within the crown has not yet gelled or set. Do not move, torque, or disturb the crown during cleanup. Following all excess removal, exposed margins may be light cured 20-40 seconds to assist restoration stabilization.

3.5 Curing, finishing and dismissal

Non-Light-Transmissible restorations (Self-Cure)

For metallic, thick or heavily opaque ceramic or composite, or restorations that otherwise impede the transmission of light, once cleanup is completed and restoration is stabilized, allow SmartCem™ 2 Cement to self-cure without disturbing for 6 minutes from start of mix. Protect restoration from contamination and movement during the setting time. Following the self-cure set, check and adjust occlusion and polish as necessary. Patient may then be dismissed.

Inadequate retention due to insufficient stabilization

- Stabilize restoration during cleanup
- Protect restoration from movement during setting
- Allow cement to self-cure without disturbing for 6 minutes

Light-Transmissible restorations (Light-Cure)

For most non-metallic, light-transmissible ceramic or composite restorations, SmartCem™ 2 Cement will be visible during cure. Once cleanup is completed and restoration is stabilized, light cure all areas of the restoration using a visible light, curing unit designed to cure CO initiated methacrylates (spectral output including 470nm), with a minimum output of 550mW/cm² for 10 seconds from each direction – buccal, lingual and occlusal. Following the light-curing, check and adjust occlusion and polish as necessary. Patient may then be dismissed.

Inadequate polymerization due to insufficient curing

- Check compatibility of curing light
- Check curing cycle
- Check curing output before each procedure

Important Technique Tips:

- When simultaneously cementing multiple single units or bridgework, it is recommended to employ the light-cure cleanup on one or two adjacent units only, allowing other units' excess to self-cure, providing ample cleanup time.
- Clean excess cement from metal instruments immediately as set cement will adhere to the instrument.

4. HYGIENE

- Cross-contamination**
 - Do not reuse single use products. Dispose in accordance with local regulations.
 - Reprocess reusable products as described below.

4.1 Cleaning

For dual-barreled syringe, remove used mixing tip and discard appropriately. Replace original syringe cap prior to storing. The dual-barreled syringe may be cleaned by scrubbing with a disposable towel soaked with hot water and soap or detergent.

4.2 Disinfection and/or sterilization

Disinfect dual-barreled syringe with a hospital-level, tuberculocidal disinfectant solution according to national/local regulations. Iodophors, sodium hypochlorite (5.25%), chlorine dioxide and dual or synergized quaternary ammoniums are approved disinfectants. Some phenolic-based agents and iodophor-based products may cause surface staining. The disinfectant manufacturer's directions should be followed properly for optimum results. Water-based disinfectant solutions are preferred.

5. LOT NUMBER AND EXPIRATION DATE

- Do not use after expiration date. ISO standard uses: "YYYY/MM."
- The following numbers should be quoted in all correspondences:
 - Reorder Number
 - Lot number
 - Expiration Date

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1.2 Composición

Urethane Dimethacrylate; Di- and Tri-Methacrylate resins; Phosphoric acid modified acrylate resin; Barium Boron FluoroAluminoSilicate Glass; Organic Peroxide Initiator; Camphorquinone (CQ) Photoinitiator; Phosphine Oxide Photoinitiator; Accelerators; Butylated Hydroxy Toluene; UV Stabilizer; Titanium Dioxide; Iron Oxide; Hydrophobic Amorphous Silicon Dioxide

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