# Zelgan®2002

## Dustfree Alginate Impression Material

DeTrey Zelgan®2002 is a dust-free alginate impression material complying with ISO 1563 and EN 21563, respectively. The powder has been especially formulated to eliminate dust during mixing.

### COMPOSITION

- Diatomaceous earth (amorphous SiO<sub>2</sub>)
- Potassium alginate
- Calciumsulfate dihydrate
- Polypropylene glycol
- Disodium orthophosphate
- Potassium hexafluorotitanate
- Sodium aluminosilicate (zeolith type A)
- Magnesium oxide
- Red iron oxide pigment
- Spearmint oil

### INDICATIONS

Dental impressions for the fabrication of

- Study models
- Models of the opposing dentition
- Working models for manufacturing of orthodontic appliances, removable prosthesis, and temporary restorations

### CONTRAINDICATIONS

Do not use in patients with a known allergy to peppermint oil or other components of the material.

### WARNINGS

Zelgan 2002 may cause allergic reactions in patients sensitised to Spearmint oil.

In isolated cases allergic reactions have been reported after impression taken with Zelgan 2002

### INTERACTIONS WITH DENTAL MATERIALS

#### Use of Tray Adhesive

For optimum adhesion, a tray adhesive for alginate materials, e.g. DeTrey® Fix, may be used.

### STEP-BY-STEP INSTRUCTIONS

#### Proportioning

Alginate impression materials are fine powders which tend to settle during transit and storage. After transferring the contents of the bag to the plastic container, always loosen the powder before use by inverting the container twice. Sufficient volume for this has been allowed for in the container. For correct proportioning, scoop powder lightly and level off; do not tap the scoop.

Sectional impression:	1 level scoop of powder	+ lower mark of water measure
Full impression:	2 level scoops of powder	+ middle mark of water measure
Large impression:	3 level scoops of powder	+ upper mark of water measure

If a thinner or thicker mix is required use slightly more or less water.

The standard powder : water ratio<sup>1</sup> of 48 ml alginate powder (2 scoops) and 34 ml water (middle mark of water measure) results in 39 ml impression material.

### Mixing

#### Machine mixing

Mix for 8 to 12 seconds (or according to manufacturer's instructions). Machine mixing will slightly shorten both working and setting times.

#### Handmixing

Mix vigorously in a suitable bowl for the time stated in the table below.

### Working and Setting Times<sup>2</sup>

Working and setting times can be modified by altering the temperature of the mixing water (see Table).

Water temperature (°C)	Mixing time (seconds)	Minimum working time (minutes)	Approx. setting time at 35 °C from start of mix (minutes)
15	30	2'05"	3'00"
23	30	1'30"	2'00"
31	30	0'50"	1'20"

The above setting times assume that the material is placed in the mouth 30 seconds from the end of the mix. In practice the time of placement may vary altering working and setting times. In general, at normal room temperature, the impression should remain in the mouth for approximately 1 minute.

### Impression Taking

Have the patient rinse the mouth with warm water. Load and seat the tray into the mouth immediately after completion of mixing and retain, under slight pressure for approximately 1 minute. Remove the impression from the mouth and rinse free from saliva and debris, under cold running water.

### Cleaning, Disinfection, Sterilisation

For disinfection immerse the impression in a 2 % aqueous solution of glutaraldehyde for 10 minutes (e. g. Impresept<sup>3</sup>). After disinfection, rinse impression under running water for approximately 15 seconds.

### Casting

It is recommended that the model be cast from the impression immediately after removal from the mouth. If this is not possible then wrap the impression in a damp napkin and keep in a sealed polyethylene bag to minimise any possible dimensional change.

Zelgan 2002 impressions will retain their accuracy for a considerable time if stored in a humidior maintaining at least 90 % relative humidity but it is preferable that the model be cast within 48 hours of taking the impression.

Compatibility with cast materials has been confirmed for Moldano blau (Heraeus Kulzer), Vel Mix Stone (Kerr), Fuji Rock EP (GC), Calspar Rapid DP, Calstone Hard DP and Calstone Super 4.

### STORAGE

Store Zelgan 2002 in a tightly closed container at temperatures between 2 and 28 °C for longer periods of time. Inadequate storage conditions will shorten the shelf life and may lead to malfunction of the product.

### BATCH NUMBER AND EXPIRY DATE

The batch number should be quoted in all correspondence which requires identification of the product.

Do not use after expiry date.

<sup>1</sup> For testing according to ISO 1563: Powder : water ratio = 14.7 g powder : 34 ml water.

<sup>2</sup> The quoted working time is based on the use of distilled water but there is no significant change if tap water of any hardness is used.

<sup>3</sup> Impresept is a trademark of 3M ESPE, Seefeld, Germany.

If you have any questions, please contact:

**Manufacturer:**  
DENTSPLY DeTrey GmbH  
De-Trey-Str. 1  
D-78467 Konstanz  
Phone (0 75 31) 5 83-0

**Distributor:**  
DENTSPLY Limited  
Hamm Moor Lane  
Addlestone, Weybridge  
Surrey KT15 2SE, UK  
Phone (0 19 32) 85 34 22

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